# AUTOCLAVE ESTERILIZADOR AUTOMATICO AUTOMATIC AUTOCLAVE STERILIZER

Manual de instrucciones

**SA-300VF**  
**SA-300VL**

## CONTENIDO

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## CONTENTS

Please read this manual carefully prior to using your SA-300VL.

Following the simple instructions contained in this manual will help ensure ease of use, trouble free operation, and a longer working life for this Autoclave.

CHAPTER 1. LOCATION OF CONTROLS

CHAPTER 2. SPECIFICATION

CHAPTER 3. PREPAEATION & INSTALLATION

CHAPTER 4. OPERATION & BASIC STEPS FOR STERILIZING

CHAPTER 5. WARNING

CHAPTER 6. APPENDIX

CHAPTER 7. SIMPLE TROUBLE SHOOTING


CHAPTER 8. MAINTENANCE REQUIREMENTS

WARRANTY

ISO-13485 CERTIFICATED FIRM

## CHAPTER 1. LOCATION

1. PRESSURE EXHAUST NEEDLE VALVE FOR LIQUID PROGRAM
2. DOOR CLOSED SCREW RING DISPLAY
3. DOOR HANDLE
4. PRESSURE / TEMPERATURE GAUGE
5. ABS COVER OF CONTROL BOX
6. POWER LAMP (INDICATION LAMP)
7. HEATING, STERILIZATION LAMP (INDICATION LAMP)
8. COMPLETE LAMP (INDICATION LAMP)
9. ACTURAL CHAMBER TEMPERATURE DISPLAY
10. SET-UP STERILIZATION TEMPERATURE
11. TEMPERATURE SET-UP FUNCTION KEY
12. STERILIZATION TIMER
13. START SWITCH
14. POWER SWITCH
15. ADJUSTABLE FOR EXHAUST

## CHAPTER 2. SPECIFICATION

- CHAMBER SIZE: 300mmØ \* 710mm (D).
- CHAMBER CAPACITY: APPROX. 50 Liter.
- OVERALL: 600 mm (W) \* 1090mm (H) \* 450 mm (D).
- POWER CONSUMPTION: 230VAC / 2.9KW / 13A.
- DESIGN TEMPERATURE: 140°C
- SAFETY DEVICE:
  - OVER HEAT (LOW WATER) CUT-OFF SWITCH
  - SAFETY VALVE AND RELEASE VALVE
  - EMERGENCY EXHAUST VALVE
  - PRESSURE PROTECTION SWITCH
  - TEMPERATURE PROTECTION SWITCH
  - TEMPERATURE CONTROL SYSTEM
- LIFE TIME: 7 YEARS

#### - STANDARD ACCESSORIES:

- Stainless Steel #304 Basket 270mmØ×308mm(D) /2pc.
- Stainless Steel #304 easy shovel / 1pc.
- 2000cc bottle / 1pc.
- Exhaust hose / 1pc.

- SAFETY APPROVAL: CE-EMC, LVD APPROVED, CE-MDD APPROVED, CHAMBER DESIGNED BY ASME / U.S.A., STANDARD

NE

### CHAPTER 3. PREPARATION and INSTALLATION

#### CHECKLIST BEFORE USE THIS AUTOCLAVE

1. Put this autoclave on the level floor, and keep the distance more than 5cm between the wall and the case. And keep this autoclave in level condition.

PS: Don't forget to press stopper on the wheel to fix the autoclave.

2. Check the elec. Power source 220v~240v.

3. You are now, ready to take the power cord out, and connect the wire to 1 Ø, 20 Amps. With ground wire and No Fuse Breaker.

PS: SA-300VL WITHOUT FUSE INSERT, SO THE NO-FUSE-BREAKER IS NECESSARY FOR SECURITY.

<<VERY IMPORTANT>>

4. Pouring about 3000cc water (DISTILLED WATER ONLY) into the chamber, (As our enclosed one and half of bottle water.)

5. Install the exhaust hose: One end to be fixed on the pipe connector "A". And another end must be fixed on the connector of the water exhaust piping system.

A. Steam exhaust hose from safety valve.


B. Pressure / water exhaust connector.

C. Power cord.

6. Push the power switch-on, if it is light on, that means the power is already stand-by now.

7. Please don't forget to check the #14 AND #1 exhaust valve (at chapter 1) must be closed, before you start the function.

8. If you can complete above preparation, She'll be more safety and she'll be your vest partner.

### CHAPTER 4. How to operation the UT150 temperature control.

Power on the machine, and wait 5 second for machine self-test.

\*Set the sterilization temperature.

1. Please press ▼ or ▲ to set sterilization temperature 121 or 135°C
2. Press SET/ENT key again

PS Please do not change following data. Because the data is for you reference. If you change the data will crush the machine. The difference machine has difference data.

1. Press SET/ENT key about 3 second. The display is PV .....A1  
SP..... 0

2. Press SET/ENT key. The display is PV ... CTL  
SP... PID
3. Press SET/ENT key. The display is PV ... AT  
SP ... oFF
4. Press SET/ENT key. The display is PV ... P  
SP ... 4.0
5. Press SET/ENT key. The display is PV ... I  
SP ... 4.
6. Press SET/ENT key. The display is PV ... D  
SP ... 10.
7. Press SET/ENT key. The display is PV ... CT  
SP ... 2.
8. Press SET/ENT key. The display is PV ... FL  
SP ... oFF
9. Press SET/ENT key. The display is PV ... BS  
SP... 1.3
10. Press SET/ENT key. The display is PV ... LoC  
SP ... 0

Press SET/ENT key about 3 second. The machine come back operation model.  
The PV value is xxx. The SP is 121 or 135 (The value is you set).

## CHAPTER 4. -1 OPERATION

\*PLEASE CHECK CHAPTER 3. AT FIRST

1. Open the door and pour water (one and half of bottle water – 3000cc) in the chamber. Then, put-in sterilized instruments in the basket.

PS: PLEASE DON'T FORGET TO PUT THE STERILIZATION INDICATORS INTO CHAMBER

2. Closed the door the turn the door handle to tight.

\*ATTENTION: DOOR HANDLE MUST BE TIGHT.

PLEASE TURN-OFF THE WATER EXHAUST VALVE AND PRESSURE RELEASE NEEDLE VALVE BEFORE YOU START THE CYCLE.


3. Push-ON the power switch, and waiting about 5 seconds until the temperature digital display in stand-by condition. Then, set-up the sterilization temperature. Please check the CHAPTER 4 for operate of temperature control.

4. Set-up the sterilization time, as you want. Please reference the chapter 6-A.

5. Press "START" switch to start program.

6. This autoclave will alarm you, when she finishes work. And please check the pressure gauge, if the chamber still in pressure condition, please:

6-1. If you want to exhaust the water in the chamber, please turn open the pressure exhaust valve.

6-2. If you want to release the chamber pressure only, please turn open the pressure release needle valve, for the liquid sterilization.

6-3-1 If you don't want to open the door soon, you can wait the chamber temperature cooling down (approx. 30 minutes).

6-3-2 If you want to open the door soon, you can turn open the pressure release needle valve 180°. The pressure will release approx. 15 minutes.


#### CAUTION

1. Please don't forget to turn-off the exhaust valve and pressure release needle valve again, when the pressure return 0 kg/cm<sup>2</sup>.

2. Please don't forget to turn-off the "START" switch to off position.

3. TO USE THE UNIT CONTINUOUSLY ARE ABSOLUTELY SURE TO OPEN THE DOOR AND ALLOW THE UNIT TO COOL OFF FOR AT LEAST 20 MINUTES.

## CHAPTER 4. -2 BASE STEP OF STERILIZATION

The door can't open, if chamber in vacuum condition.

1. Please turn pressure exhaust needle valve to release chamber vacuum.

2. after the chamber vacuum was released, then the door can be opened.

Caution:

1. Check the pressure gauge return to "zero" position, before the door opened.

2. If you need do the next run. please turn off/on the power switch to reset system again.

3. To use the unit continuously are absolutely sure to open the door and allow the unit to cool off for at least 20 minutes.

## CHAPTER 5. WARNING

1. Please check the pressure gauge, when you open the door. (PLEASE DON'T OPEN THE DOOR, IF THE PRESSURE OVER THAN 0 kg/cm<sup>2</sup>).

2. Please use the single power source, at the same time. And please don't forget to device a 1Ø 20Amps. NO FUSE BREAKER between the autoclave and elec. Power source. IT IS VERY IMPORTANT FOR THE SAFETY REQUEST.

3. "DISTILLED WATER" is necessary for SA-300VL.
4. The "LOW WATER" a long audio sound that is means the water in chamber is not enough.
- 4-1. Please check the exhaust valve was closed complete or not.
- 4-2. If you set-up the sterilization time over than 40 minutes, please pour 3000cc water (1.5 bottle) or more.
5. The door holder must be tight, when you close the door.
6. Please keep the chamber and gasket clean.
7. Please don't forget to put in the sterilization indicators. And don't forget to check the sterilization indicators, after sterilization.
8. Please attend the high temperature on the cover of the door. When she is working.
9. Every week, must clean the chamber and basket. To keep the autoclave in good condition.
10. The exhaust valve also can be EMERGENCY PRESSURE RELEASE device. When the autoclave out of control, or the user want to stop the working program, immediately. You can turn-open this valve about 50°-90°. Until the pressure return to 0 kg/cm<sup>2</sup>, you can open the door.
11. THE SIGN OF CAUTION AND INDICATION:

CAUTION: Please check our Operation manual.

CAUTION: Hot, Please don't Touch it.

INDICATION: Electricity earth Protection.

INDICATION POWER – ON

INDICATION POWER - OFF

PS: SOME TIME, AFTER THE CHAMBER COOLING RETURN (AFTER WORKING), THE CHAMBER WILL HAVE A VACUUM CONDITION. PLEASE TURN-OPEN THE EXHAUST VALVE. THEN, THE VACUUM CONDITION WOULD BE RELEASED IMMEDIATELY.

12. Anytime when you push-on the power switch, please wait about 5 seconds. And then you can set-up the sterilization temperature and sterilization time. Because the temperature control switch need 5 seconds for auto-research all the function and programs.
13. If you push the "SET" bottom over than 3 seconds, the display will be jumped to other functions. Please push this button 3 seconds again. And the display will come back to normal condition.

## CHAPTER 6. APPENDIX- A

FOR REFERENCE OF PROGRAM SETTING-UP:

STERILIZATION TEMPERATURE SELECTION RANGE: 105°C to 136°C

STERILIZATION TIME SELECTION RANGE: 0 to 60 MINUTES.

#121°C Approx. 1.1 kg/cm<sup>2</sup>

#134°C Approx. 2.2 kg/cm<sup>2</sup>

# WRAPPED STERILIZATION TIME      15minutes for 134°C

30minutes for 121°C

# UNWRAPPED STERILIZATION TIME      04minutes for 134°C

22minutes for 121°C

# LIQUID STERILIZATION TIME      40minutes for 121°C

#FOR REFERENCE: 1KG/CM<sup>2</sup> = 0.98 bar = 14.2 psi.

STORAGE ENVIRONMENT: TEMPERATURE: -10°C~+50°C / HUMIDITY: ≤80%

WORKING ENVIRONMENT: TEMPERATURE: 1°C~+40°C / HUMIDITY: ≤80%

TRANSPORTATION ENVIRONMENT: TEMPERATURE: -10°C~+60°C / HUMIDITY: ≤80%

## CHAPTER 6 APPENDIX –B

\*RE-PACKING INSTRUCTIONS: TO PACK AUTOCLAVE INTO CARTON (PLEASE REFER CHAPTER 1.LOCATION.)

1. Drain all water from the autoclave.
2. Remove power cord and place all items into the chamber. Close door secure filmy.
3. Cover the autoclave with the plastic bag.
4. Place the bottom of the wooden-case onto the floor.
5. Move and place the autoclave on bottom piece. And ensuring that these four castors are fixed already.
6. Place wall pieces of the wooden case to connection the bottom piece of carton, one by one.
7. Place the heavy wooden bar to fix the top of the autoclave, between each wall pieces.
8. Putting on the top cover of the wooden –case.
9. Please don't forget to write down the shipping mark and products name, model No. on the same time.

IMPORTANT: We strongly recommend that all packing material including the carton be kept for re-use in the event that the unit needs to be transported at ant time.

## CHAPTER 7. TROUBLE SHOOTING

POWER LAMP IS NOT LIGHT ON:

1. Power supply is not properly connected.  
/ 1. Try to connect power supply until power indicator goes on.
2. No Fuse breaker broken.  
/ 2. Replace the No Fuse Break. (Please check CHAPTER 1.)
3. Bulb or switch broken. / 3. Replace bulb or switch.

TWO INDICATION LAMPS LIGHT-ON AND ALARM:

1. Water is not enough. / 1. Pull the water in the chamber. (Please pour water in, with the enclosed 3000cc bottle – one and half of bottle water) or 1-1. check the exhaust valve is not be closed complete, or the valve broken.

PS: If you set-up the sterilization time over than 40 minutes, please pour 3000cc water (one and half of bottle water)

PRESSURE CAN'T UP:

1. Steam solenoid unclean. / 1. Replace clean the filter.
2. Heater broken / 2. Replace the heater
3. Temperature protection switch broken / 3. Replace the switch.
4. Pressure gauge broken / 4. Replace the gauge.
5. Needle valve open / 5. Close needle valve.

THE STERILIZATION INDICATOR COLOR UNCHANGED:

1. The temperature can't coming-up / 1.Please follow above mention.
2. The steam trap broken / 2. Replace the steam trap.

THE INDICATION LAMPS CAN'T LIGHT-ON:

1. If the function cycle is order / 1. Replace the bulbs

PS: ANY PROBLEMS, PLEASE CHECK THE WIRE CONNECTION AT FIRST.

## CHAPTER 8. MAINTENANCE REQUIREMENTS

**DAILY:**

WIPE THE INSIDE OF THE CHAMBER AND THE INSIDE OF THE DOOR WITH NON-LINT CLOTHE SUCH AS A WETTEX.

**WEEKLY:**

- CLEAN THE STERILIZATION BASKET, AND DOWN TRAY.

- CLEAN THE FILTER IN THE CHAMBER.
- CLEAN THE CHAMBER.

#### MONTHLY:

- CHECK THE ELECTRICITY POWER CORD, BREAKER AND PLUG SURFACE IN GOOD CONDITION. IF NOT, PLEASE CALL THE ENGINEER.
- PUSH ON THE POWER SWITCH, AND CHECK THE FAN AT REAR SIDE IS WORKING. IF NOT, PLEASE CALL THE ENGINEER.
- CLEAN THE CHAMBER AND PIPING SYSTEM WITH "CHAM-MATE".

#### YEARLY: (CHECK BY ENGINEER)

- CHECK THE WELDING LINE OF THE CHAMBER.
- CHECK THE DOOR LOCK SCREW.  
DOOR HOLDER SCREW AND HOLDER BODY.
- CHECK THE ELECTRIC WIRE SYSTEM, BREAKER, AND CONNECTORS.
- CHECK THE TUBING SYSTEM.
- CHECK THE SAFETY VALVES, PRESSURE SWITCH, AND LEAKAGE TEST.
- CHECK THE INDICATION LAMPS, AND ALL CYCLE FUNCTION.
- CHECK THE DOOR GASKET. (SUGGESTION: REPLACE IT PER YEAR.)

#### WARRANTY

Quirumed offers the guarantee provided by law against all manufacturing defects from the date of purchase. You should therefore retain the invoice or purchase.

When using non-original parts or modifies the product immediately void.

Misuse or wear of product is not considered a manufacturing defect and hence not included in the warranty.

Nor a bad unpacking the product included such cuts made with cutter or casualties.

For any matter relating to the warranty please contact your dealer.

Distribuido por: Quirumed, S.L. - C.Moroder, 3 Pol. Ind. III Moncada 46113 - Valencia (Spain)  
Tel.: 96 338 38 33 - Fax: 96 338 38 34 - [www.quirumed.com](http://www.quirumed.com)

APEX MEDICAL S.L. - Elcano 9, 6a planta 48008 Bilbao, Vizcaya SPAIN

Manufacturer: STURDY INDUSTRIAL CO.,LTD. (ISO 13485 CERTIFICATE FIRM.)  
Sturdy Industrial Co., Ltd. - No. 168, Sec. 1, Zhongxing Road, Wugu Dist., New Taipei City 24872, Taiwan