# AUTOCLAVE ESTERILIZADOR AUTOMATICO AUTOMATIC AUTOCLAVE STERILIZER

Manual de instrucciones

**SA-232V**

## CONTENIDO

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Distribuidor/*Distributor*: \_\_\_\_\_

## CONTENTS

Please read this manual carefully prior to using your new Autoclave. Following the simple instructions contained in this manual will help ensure ease of use, trouble free operation, and a longer working life for your autoclave.

CHAPTER 1. LOCATION .....	8
CHAPTER 2. SPECIFICATION .....	8
CHAPTER 3. PREPARATION.....	9
CHAPTER 4. OPERATION .....	9
CHAPTER 4-1. BASIC STEP OF STERILIZATION.....	10
CHAPTER 5. WARNING.....	10
CHAPTER 6. APPENDIX -A .....	11
CHAPTER 6. APPENDIX -B .....	11
CHAPTER 7. TROUBLE SHOOTING .....	11
CHAPTER 8. MAINTENANCE REQUIREMENTS .....	12

### CHAPTER 1. LOCATION (page 2)

1. Safety valve
2. Door handle
3. Door closed screw
4. Pressure gauge
5. Chamber pressure exhaust valve (Exhaust flowing speed adjusted by knob)
6. Pressure control switch (Adjustable)
7. Power indicator
8. Heating-up indicator
9. Sterilization indicator
10. Complete indicator
11. Program start switch
12. Sterilization time: 4 minutes
13. Sterilization time: 15 minutes
14. Sterilization time: 30 minutes
15. Sterilization time selection switch
16. Main power switch

### CHAPTER 2. SPECIFICATION

CHAMBER SIZE	SA-232V 230mmØ X 400mm(DEPTH)
OVERALL SIZE	SA-232V 380mm(W) x 600mm(H) x 365mm(D)
POWER CONSUMPTION	SA-232V 230VAC /1.2 KW/ / 5.3A
SAFETY DEVICE:	
- PRESSURE SAFETY VALVE	
- LOW WATER POWER CUT-OFF	
- TEMPERATURE PROTECT SWITCH	
- PRESSURE CONTROL SWITCH	
- EMERGENCY EXHAUST KNOB	
- DOOR SAFETY CATCHER	
DESIGN TEMPERATURE	140°C
STANDARD ACCESSORIES	STERILIZATION BASKET x 1. (st. #304) HEATER COVER x 1 (st. #304)
LIFE TIME	7 YEARS

## CHAPTER 3. PREPARATION

### CHECKLIST FOR UNPACKING INSTRUCTIONS:

TO UNPACK AUTOCLAVE FROM CARTON

1. Remove straps.
2. Remove top of carton.
3. Remove center of piece (walls) of the carton.
4. Remove foam-packing pieces.

### CHECKLIST BEFORE USE THIS AUTOCLAVE

**NOTE:** PLEASE CHECK THE ACCESSORIES ENCLOSED, AT FIRST.

1. Put this autoclave on the table, and keep the distance more than 5cm between the wall and the case. And keep this autoclave in level condition.
2. Check the elec. power source must be 230VAC.  
230V AC power cord type: 0.75mm<sup>2</sup> x 3C TUV
3. Pouring the water (DISTILLED WATER or FILTER WATER ONLY) into the chamber, until the water level FULL COVER the heater (about 1000 cc for SA-232V).

**CAUTION:** PLEASE DON'T OVER THAN THE HEATER COVER.

4. You are now, ready to take the power cord out from autoclave, and plug the female end into the standard grounded AC outlet.

**NOTE:** POWER SOCKET 230VAC /10A

**CAUTION:** Before plug the power cord, please turn-off the power switch.

5. Push-on the power switch; check the power lamp, if it is light on, that means the power is already stand-by now.

## CHAPTER 4. OPERATION

PLEASE CHECK CHAPTER 3 AT FIRST

1. Open the door and pouring water (SA-232V approx. 1000cc) into the chamber. Then, put the sterilized instruments in the sterilization basket.

**NOTE:** PLEASE DON'T FORGET TO PUT THE STERILIZATION INDICATORS INTO CHAMBER

2. Closed the door and turn the handle screw to tight.

**ATTENTION:** DOOR HANDLE MUST BE TIGHT.

**PLEASE TURN-CLOSED THE WATER EXHAUST VALVE COMPLETELY**

3. Push on the power switch (The power indicator light on).
4. Set-up the pressure control knob 0.9 kgf/cm<sup>2</sup> (115°C) to 2.1 kgf/cm<sup>2</sup> (134°C).
5. Set-up the sterilization time:

Set up the sterilization time. Please note: The sterilization time was set up in the programs 4 minutes / 15 minutes / 30 minutes. So the user can select the sterilization time by selection switch. And then, push on the "START" switch on the same time; the heating lamp will light on.

AS TO TIME SET, PLEASE REFER THE ENCLOSED APPENDIX.

6. After the chamber pressure up to the pressure you set up, the sterilization indicator will light on. Meanwhile the sterilization timer starts work when the sterilization indicator light off (and the "COMPLETE" indicator light on) that means the sterilization cycle is completed.
7. 40 seconds buzzer will report you, "MY DEAR MASTER! YOUR COMMAND HAVE BEEN COMPLETED".
8. If you want to open the door immediately, after cycle completed. Please turn open the exhaust valve to 25° for liquid sterilization, 50° for normal sterilization.

### CAUTIONS

1. Please don't forget to turn-off the exhaust valve when the pressure returns 0 kgf/cm<sup>2</sup>
2. Before opening the door ensure the pressure gauge is at "zero" (0) position.
3. Please don't forget check the sterilization indicator after sterilization cycle.

4. To use the unit continuously are absolutely sure to open the door and allow the unit to cool off for at least 20 minutes.

## CHAPTER 4-1. BASIC STEP OF STERILIZATION

TURN CLOSED THE EXHAUST VALVE (MUST CLOSED COMPLETELY).

POUR WATER INTO CHAMBER; THE WATER LEVEL COVERS THE HEATER COMPLETELY.

PUT THE STERILIZATION INSTRUMENTS INTO CHAMBER

1. SET UP STERILIZATION PRESSURE
2. SELECT STERILIZATION TIME
3. CLOSED DOOR.

PUSH START SWITCH (START WORKING) THEN MAIN HEATER ON STERILIZATION TIMER START WORKING, WHEN PRESSURE ARRIVE YOU SET-UP

BUZZER ON, WHEN STERILIZATION TIMER IS TIME OUT

TURN THE ADJUSTABLE EXHAUST SWITCH FOR EXHAUST WATER OF CHAMBER

**CAUTION:** CHECK THE PRESSURE GAUGE RETURN TO "ZERO" POSITION, BEFORE THE DOOR OPENED.

## CHAPTER 5. WARNING

1. Anytime check the pressure gauge, if the pressure over than Okgf/cm2. Please don't open door.
2. "USE ONLY DISTILLED WATER" or not, we can't offer the quality guarantee service.
3. Please pour the water into the chamber before sterilization. Or not, the low-water protector will cut-off the elec. power. Please TURN-OFF the power switch, and pour water into chamber. Then, close the door and TURN-ON the power switch again.
4. Please attend the high temperature on the top of the case. When she is working.
5. Ensure the heater element inside the chamber is covered with water prior to commencement of the sterilization cycle.
6. Use only a dedicated power supply.
7. The silicon rubber gasket on the inner door and the front of the chamber should be kept clean.
8. The filter valve inside the chamber should be cleaned at least every season.
9. We recommend use of chemical indicator strips as a check for sterility. These strips may also be kept as a record of sterilization.
10. Do not over load or overfill the water into the water reservoir.
11.  This label is meaning this place is hot will serious injury of the operator or other persons.
12.  This label is meaning that this place is ground.
13.  This label is meaning that indicator for power off.
14.  This label is meaning that indicator for power on.
15.  Please careful read the operation manual.
16.  The parts need to recycle cannot throw away.
17.  2008-01 Date of Manufacture. EX: 2008-01 is mean January / 2008.
18.  Authorized Representative at europe

19.  Manufacturer.

20. Storage environment: Temperature: -10°C ~+50°C / Humidity::80% RH

21. Working environment: Temperature: 1°C ~+40°C / Humidity::80% RH

22. Transportation environment: Temperature: -10°C ~+60°C / Humidity≤80% RH

## CHAPTER 6. APPENDIX -A

### SETTING UP THE STERILIZATION TIME:

**PRESSURE ON 1.2 kgf/cm<sup>2</sup> (APPROX. 121°C / 250°F)**

UNWRAPPED Set-up 15 minutes sterilization time.

WRAPPED Set-up 30 minutes sterilization time

LIQUID Set-up 30 minutes sterilization time.

For liquids sterilization : Please turn adjusted the exhaust valve on 20°-30°.

For normal sterilization : Please turn adjusted the exhaust valve on 50°.

- FOR REFERENCE: 1 kgf/cm<sup>2</sup> = 0.98 bar = 14.2 psi.

- Sterilization instruments must be made by metal without plate, or autoclavable plastic or silicon rubber material.

- Liquid sterilization is available, and SA-202V and SA-232V is special design for laboratory and liquid sterilization.

## CHAPTER 6. APPENDIX -B

RE-PACKING INSTRUCTIONS: TO PACK AUTOCLAVE INTO CARTON (PLEASE REFER CHAPTER 1. LOCATION.)

1. Drain all water from the autoclave.

2. Remove power cord and place all items into the chamber. Close door secure firmly.

3. Cover the autoclave with the plastic bag.

4. Place the bottom of the carton containing the packing piece NO.1 onto the floor.

5. Place the autoclave into packing piece NO.1 ensuring that the door of the autoclave is in the position marked "door".

6. Place center piece (wall) of carton into the bottom of the carton.

7. Place packing piece NO.2 in front of door.

8. Place packing piece Nos. 3, 4, and 5 (all the same size) along both sides and rear of autoclave.

9. Place packing piece No.6 on top of autoclave.

10. Place top of carton over center piece (wall) of carton.

11. The carton together with heavy twine or a similar strong material.

**IMPORTANT:** We strongly recommend that all packing material including the carton is kept for re-use in the event that the unit needs to be transported at any time.

## CHAPTER 7. TROUBLE SHOOTING

### • POWER LAMP IS NOT LIGHT ON:

1. Power supply is not properly connected.

Solution: Try to connect power supply until power indicator goes on.

2. Fuse broken

Solution: Replace fuse or replace fuse socket.

3. LED broken.

Solution: Replace the PCB or LED.

### • STERILIZATION LAMP IS NOT LIGHT ON:

1. LED Broken.

Solution: Replace the LED

2. Relay broken

Solution: Replay relay (on PCB)

3. Pressure switch broken

Solution: Replace the pressure switch.

• **LOW WATER AND ELEC. POWER CUT OFF:**

1. The exhaust valve don't closed complete

Solution: Turn close complete.

2. For pouring water into the chamber or the water not enough.

Solution: Pull the water into the chamber; the water level must cover the heater (Don't over than the base tray).

3. Thermo switch broken

Solution: Replace the thermo switch.

• **PRESSURE CAN'T UP:**

1. The exhaust valve don't closed complete

Solution: Turn close complete

2. Heater broken.

Solution: Replace main heater

**Note:** ANY PROBLEMS PLEASE CHECK THE WIRE CONNECTION AT FIRST.

## CHAPTER 8. MAINTENANCE REQUIREMENTS

### DAILY:

- Wipe the inside of the chamber and the inside of the door with non-lint clothe.

- Ensure the exhaust hose connection fixed on the autoclave and other end.

### WEEKLY:

- Clean the basket and heater cover.

- BIO-test the sterilization result. If the result is fail, Please call the engineer to check it.

### MONTHLY:

- Check the power cord and plug. If the plug and socket have been melted. Please call engineer to check the power source and the electric wire system.

- Clean the chamber and piping system with piping cleaner.

- Pull three times for safety valve.

### YEARLY: (Check by engineer)

- Check the electric wire system, fuse, and connectors.

- Check the tubing system.

- Check the indication lamps, and all cycle function.

- Check the door gasket. (Suggestion replace it per year.)

## WARRANTY (page 7)

This product are in guaranteed to be from the defects in materials and workmanship for one (1) year under normal use from the date of purpose.

This warranty does not apply to any product damaged by accident, misuse abuse, neglect, improper line voltage, drop, fire, flood or if the products were altered or repaired by anyone other than qualified service personnel.

The liability of company is limited to repair or replacement and under no circumstances shall be liable for any collateral consequential damages or loss. This guarantee specifically excludes the expendables and consumable.

All warranty claims must be directed to the distributors or agents authorized by company responsible for the sales of this instrument. The customers are responsible for shipping expense.

Distribuido por: Quirumed, S.L. - C.Moroder, 3 Pol. Ind. III Moncada 46113 - Valencia (Spain)  
Tel.: 96 338 38 33 - Fax: 96 338 38 34 - [www.quirumed.com](http://www.quirumed.com)

APEX MEDICAL S.L. - Elcano 9, 6a planta 48008 Bilbao. Vizcaya SPAIN

Manufacturer:  
Sturdy Industrial Co., Ltd. - No. 168, Sec. 1, Zhongxing Rd., Wugu District, New Taipei City 24872, Taiwan