

No.		Item	Parameter for Prime 30L, 56L			Parameter			Parameter			Parameter			Remark / Reference		
			Prime 30L Economic	Prime 30L Standard	Prime 30L Advance	Prime 56L Economic	Prime 56L Standard	Prime 56L Advance									
1			Pressure Vessel characteristics														
1.01	1	Chamber Capacity	30 Liters / 7.9 Gallons	30 Liters / 7.9 Gallons	30 Liters / 7.9 Gallons	56 Liters / 14.8 Gallons	56 Liters / 14.8 Gallons	56 Liters / 14.8 Gallons				Prime 30 PED data(PED TCF C301-1.1)			Prime 56 PED data(PED TCF C380-1.1)		
1.02	2	Chamber Size	Ø 300 x 452 (mm) / Ø 11.8 x 17.8 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Ø 300 x 452 (mm) / Ø 11.8 x 17.8 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Ø 380 x 452 (mm) / Ø 11.8 x 17.8 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Ø 380 x 522 (mm) / Ø 15 x 20.6 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Ø 380 x 522 (mm) / Ø 15 x 20.6 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Ø 380 x 522 (mm) / Ø 15 x 20.6 (inch) Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)				Prime 30 PED data(PED TCF C301-1.1)			Prime 56 PED data(PED TCF C380-1.1)		
1.03	3	Pressure Vessel Material	Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Chamber : SUS 316L Stainless Steel , Thinness 2mm Door : SUS 316 Stainless Steel , Thinness 2.5mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)	Chamber : SUS 316L Stainless Steel , Thinness 2.5mm Door : SUS 316 Stainless Steel , Thinness 3mm Roughness of surface : Polish Gas tungsten arc welding(GTAW)				Prime 30 Chamber drawing			Prime 56 Chamber drawing		
1.04	4	Design Pressure	3.0 bar	3.0 bar	3.0 bar	3.0 bar	3.0 bar	3.0 bar				PED Car.					
1.05	5	Design Temperature	144℃	144℃	144℃	144℃	144℃	144℃				PED Car.					
1.06	6	Working Pressure / Temperature	(-1.0)~ 2.1 bar /0~135℃	(-1.0)~ 2.1 bar /0~135℃	(-1.0)~ 2.1 bar /0~135℃	(-1.0)~ 2.1 bar /0~135℃	(-1.0)~ 2.1 bar /0~135℃	(-1.0)~ 2.1 bar /0~135℃				Prime 30/56 specification(DHF 4-2-1)					
1.07	7	Chamber Design Tec.	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS	Chamber Type : Single chamber w/o jacket Resistant to oxidation, high temperature and high pressure. Gas tungsten arc welding(GTAW) RoHS				Prime 30 Chamber drawing			Prime 56 Chamber drawing		
1.08	8	Chamber Door Tec.	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation	Chamber Door structure : Quick locking cover Concave type Resistant to oxidation				Prime 30 PED data(PED TCF C301-1.1)			Prime 56 PED data(PED TCF C380-1.1)		
1.09	9	Insulation system	Material : Melamine Foam + AL tape + iron tape thickness 14 mm	Material : Melamine Foam + AL tape + iron tape thickness 14 mm	Material : Melamine Foam + AL tape + iron tape thickness 14 mm	Material : Melamine Foam + AL tape + iron tape thickness 14 mm	Material : Melamine Foam + AL tape + iron tape thickness 14 mm	Material : Melamine Foam + AL tape + iron tape thickness 14 mm				Prime 30/56 BOM list(DHF 4-2-5.a)					
1.10	10	Design life	7 years	7 years	7 years	7 years	7 years	7 years				Prime 30/56 Life test(DHF 4-3-5)					
2			Overall Characteristics														
2.01	11	Overall Dimension	584 (W) x 480 (H) x 670 (D) (mm) 23 (W) x 18.9 (H) x 26.4 (D) (inch)	584 (W) x 480 (H) x 670 (D) (mm) 23 (W) x 18.9 (H) x 26.4 (D) (inch)	584 (W) x 480 (H) x 670 (D) (mm) 23 (W) x 18.9 (H) x 26.4 (D) (inch)	720(W) x 535(H) x 710(D) (mm)	720(W) x 535(H) x 710(D) (mm)	720(W) x 535(H) x 710(D) (mm)				Prime 30/56 Installation plan(DHF 4-2-2.c.i)					
2.02	12	Carcass Material	Powder coating Thickness 1mm	Powder coating Thickness 1mm	Powder coating Thickness 1mm	Powder coating Thickness 1mm	Powder coating Thickness 1mm	Powder coating Thickness 1mm				Prime 30 Carcass drawing			Prime 56 Carcass drawing		
2.03	13	Preset Program	15 Programs	17 Programs	17 Programs	17 Programs	17 Programs	17 Programs				Prime 30/56 specification(DHF 4-2-1)					
2.04	14	Sterilization Temperature	105~135℃, Adjustable	105~135℃, Adjustable	105~135℃, Adjustable	105~135℃, Adjustable	105~135℃, Adjustable	105~135℃, Adjustable				Prime 30/56 specification(DHF 4-2-1)					
2.05	15	Temperature Sensor	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure	PT100 x 2 pcs, chamber · folting K type x 3 pcs, steam generator · B-heater · exhaust filter (optional) HMI Digital show temperature, pressure				Prime 30/56 specification(DHF 3-1-1)					
2.06	16	Display	Pressure gauge grossce indicator (green, blue, blue flashing, red)	Pressure gauge grossce indicator (green, blue, blue flashing, red)	Pressure gauge grossce indicator (green, blue, blue flashing, red)	Pressure gauge grossce indicator (green, blue, blue flashing, red)	Pressure gauge grossce indicator (green, blue, blue flashing, red)	Pressure gauge grossce indicator (green, blue, blue flashing, red)				Prime 30/56 specification(DHF 3-1-1)					
2.07	17	Operation Interface	HMI , full color 5" touch panel (screen) Under 3000M	HMI , full color 5" touch panel (screen) Under 3000M	HMI , full color 5" touch panel (screen) Under 3000M	HMI , full color 5" touch panel (screen) Under 3000M	HMI , full color 5" touch panel (screen) Under 3000M	HMI , full color 5" touch panel (screen) Under 3000M				Prime 30/56 specification(DHF 4-2-1)					
2.08	18	Altitude Usage	Under 3000M	Under 3000M	Under 3000M	Under 3000M	Under 3000M	Under 3000M				Prime 30/56 specification(DHF 4-2-1)					
2.09	19	Net Weight	68 kg	73.5 kg	75 kg	100 kg	105.5 kg	107 kg				Prime 30/56 Packing data(DHF 4-2-3)					
2.10	20	Power Supply	230V ± 10%	230V ± 10%	230V ± 10%	230V ± 10%	230V ± 10%	230V ± 10%				Prime 30/56 specification(DHF 4-2-1)					
2.11	21	Current	12.6 A (ampere)	13.1 A (ampere)	13.1 A (ampere)	12.6 A (ampere)	13.1 A (ampere)	13.1 A (ampere)				Prime 30/56 specification(DHF 4-2-1)					
2.12	22	Fuses	15 A	15 A	15 A	15 A	15 A	15 A				Prime 30/56 specification(DHF 4-2-1)					
2.13	23	power code	quick socket 16A, European style / USA type	quick socket 16A, European style / USA type	quick socket 16A, European style / USA type	quick socket 16A, European style / USA type	quick socket 16A, European style / USA type	quick socket 16A, European style / USA type				Prime 30/56 specification(DHF 4-2-1)					
2.14	24	Power Consumption	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W	Main heater 2000W for steam generator Bend heater 900W Pipe heater 43W				Prime 30/56 specification(DHF 4-2-1)					
2.15	25	Pre-Vacuum	No	Yes	Yes	No	Yes	Yes				Prime 30/56 specification(DHF 4-2-1)					
2.16	26	Post-Vacuum	No	Yes	Yes	No	Yes	Yes				Prime 30/56 specification(DHF 4-2-1)					
2.17	27	Pressure & Temperature Control	Yes	Yes	Yes	Yes	Yes	Yes				Prime 30/56 specification(DHF 4-2-1)					
2.18	28	Filter	air filter for lifted vacuum	air filter for lifted vacuum	air filter for lifted vacuum	air filter for lifted vacuum	air filter for lifted vacuum	air filter for lifted vacuum				Prime 30/56 specification(DHF 4-2-1)					
2.19	29	Digital Reorder	Built in memory 2 GB for *** cycle data USB max 64GB	Built in memory 2 GB for *** cycle data Thermal Printer USB max 64GB	Built in memory 2 GB for *** cycle data Thermal Printer USB max 64GB	Built in memory 2 GB for *** cycle data USB max 64GB	Built in memory 2 GB for *** cycle data Thermal Printer USB max 64GB	Built in memory 2 GB for *** cycle data Thermal Printer USB max 64GB				Prime 30/56 specification(DHF 4-2-1)					
2.20	30	Connect type	USB to wifi USB to scanner USB to flash memory	RJ-45 for internet USB to wifi USB to scanner USB to flash memory	RJ-45 for internet USB to wifi USB to scanner USB to flash memory	USB to wifi USB to scanner USB to flash memory	RJ-45 for internet USB to wifi USB to scanner USB to flash memory	RJ-45 for internet USB to wifi USB to scanner USB to flash memory				Prime 30/56 specification(DHF 4-2-1)					
2.21	31	water trunk	4.2 liter for recycle use	4.2 liter for recycle use	4.2 liter for recycle use	16.1 liter for recycle use	16.1 liter for recycle use	16.1 liter for recycle use				Prime 30/56 specification(DHF 3-1-1)					
2.22	32	Protection	Overheat Overpressure power Response System Four Door Locking System	Overheat Overpressure power Response System Four Door Locking System Free from Trouble come Anar Chaining	Overheat Overpressure power Response System Four Door Locking System Free from Trouble come Anar Chaining	Overheat Overpressure power Response System Four Door Locking System	Overheat Overpressure power Response System Four Door Locking System Free from Trouble come Anar Chaining	Overheat Overpressure power Response System Four Door Locking System Free from Trouble come Anar Chaining				Prime 30/56 specification(DHF 4-2-1)					
2.23	33	Condenser	NA	NA	NA	NA	NA	NA									
2.24	34	Type of piping material of the steam and water system	Stainless Steel 304+ Copper+ silicon+PTFE	Stainless Steel 304+ Copper+ silicon+PTFE	Stainless Steel 304+ Copper+ silicon+PTFE	Stainless Steel 304+ Copper+ silicon+PTFE	Stainless Steel 304+ Copper+ silicon+PTFE	Stainless Steel 304+ Copper+ silicon+PTFE				Prime 30/56 specification(DHF 3-1-1)					
3			CONTROL SYSTEM for hardware														
3.01	35	Control Method	Microcomputer + HMI	Microcomputer + HMI	Microcomputer + HMI	Microcomputer + HMI	Microcomputer + HMI	Microcomputer + HMI				Prime 30/56 specification(DHF 4-2-1)					
3.02	36	Control Panel	5 inch / Touch panel	5 inch / Touch panel	5 inch / Touch panel	5 inch / Touch panel	5 inch / Touch panel	5 inch / Touch panel				Prime 30/56 specification(DHF 4-2-1)					
3.03	37	connect prot	USB x 3 pcs	USB x 4 pcs RJ-45 for internet	USB x 4 pcs RJ-45 for internet	USB x 3 pcs	USB x 4 pcs RJ-45 for internet	USB x 4 pcs RJ-45 for internet				Prime 30/56 specification(DHF 4-2-1)					
3.04	38	Loading Capacity										Prime 30/56 Operation Manual(DHF 4-2-2.a)					
4			CONTROL SYSTEM for Software														
4.01	39	Language	English	English	English	English	English	English				Prime 30/56 specification(DHF 3-1-1)					
4.02	40	Minimum display parameters	Temp:0.1℃ Pres:0.001bar	Temp:0.1℃ Pres:0.001bar	Temp:0.1℃ Pres:0.001bar	Temp:0.1℃ Pres:0.001bar	Temp:0.1℃ Pres:0.001bar	Temp:0.1℃ Pres:0.001bar				Prime 30/56 specification(DHF 3-1-1)					
4.03	41	Program	15 Programs	17 Programs	17 Programs	15 Programs	17 Programs	17 Programs				Prime 30/56 specification(DHF 4-2-1)					

Parameter for Prime 30L, S6L																	
No.		Item	Parameter				Parameter				Parameter				Remark / Reference		
			Prime 30L Economic		Prime 30L Standard		Prime 30L Advance		Prime S6L Economic		Prime S6L Standard		Prime S6L Advance				
4.04	42	sterilization Programs	1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mixed 13. Latex 14. Flash				Prime 30/S6 specification(DHF 4-2-1)
4.05	43	test programs	NA		15. Bowie & Dick Test 16. Leakage		15. Bowie & Dick Test 16. Leakage		NA		15. Bowie & Dick Test 16. Leakage		15. Bowie & Dick Test 16. Leakage				Prime 30/S6 specification(DHF 4-2-1)
4.06	44	other program	17. Cleaning		17. Cleaning		17. Cleaning		17. Cleaning		17. Cleaning		17. Cleaning				Prime 30/S6 specification(DHF 4-2-1)
4.07	45	Program function	Calculation of F0 value countdown curve chart		Calculation of F0 value countdown curve chart		Calculation of F0 value countdown curve chart		Calculation of F0 value countdown curve chart		Calculation of F0 value countdown curve chart		Calculation of F0 value countdown curve chart				Prime 30/S6 specification(DHF 4-2-1)
4.08	46	Functions	ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging		ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging		ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging		ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging		ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging		ID management ( Adm 1, super user 4, nomor user 15) Calculation of F0 value Calibration ( temperature · pressure · Altitude) Reservation Start Service Remind Day and Time Setting Error code message indication Loads Tracking Management Free from Troublesome Agar Clogging				規格書再修改，沒有提到這個價位內容
4.09	47	Cycle Time	1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning		1. Liquid 1 2. Liquid 2 3. Solid 1 4. Solid 2 5. Agar 6. Dissolution 7. user 1 8. User 2 9. Dry only 10. Waste Liquid 11. Waste Solid 12. Waste Mix 13. Latex 14. Flash 15. Bowie & Dick Test 16. Leakage 17. Cleaning				Cycle Time
5		Door Lock SYSTEM															
5.01	48	Pressure lock	automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"		automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"		automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"		automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"		automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"		automatic locking control pressure -0.09 bar ~ +0.09 bar Four Door Locking System "door switch, lock switch, solenoid lock"				Prime 30/S6 specification(DHF 3-1-1)
5.02	49	Automatic Lock	solenoid lock		solenoid lock		solenoid lock		solenoid lock		solenoid lock		solenoid lock				Prime 30/S6 BOM list(DHF 4-2-5.a)
6		HEATING SYSTEM															
6.01	50	Sterilization heating	Steam generator		Steam generator		Steam generator		Steam generator		Steam generator		Steam generator				Prime 30/S6 specification(DHF 3-1-1)
6.02	51	Power	2900 w		2900 w		2900 w		2900 w		2900 w		2900 w				Prime 30/S6 specification(DHF 4-2-1)
6.03	52	Bend heater	900 w		900 w		900 w		900 w		900 w		900 w				Prime 30/S6 specification(DHF 4-2-1)
7		Pipe and water system															
7.01	53	Type of piping material of the steam	Stainless Steel 304+ Copper+PTFE		Stainless Steel 304+ Copper+PTFE		Stainless Steel 304+ Copper+PTFE		Stainless Steel 304+ Copper+PTFE		Stainless Steel 304+ Copper+PTFE		Stainless Steel 304+ Copper+PTFE				Prime 30/S6 specification(DHF 3-1-1)
7.04	56	water trunk	4.2 liter, PP material water level high and low of clean water tank		4.2 liter, PP material water level high and low of clean water tank		4.2 liter, PP material water level high and low of clean water tank		16.1 liter,Stainless Steel 304 material water level high and low of clean water tank		16.1 liter,Stainless Steel 304 material water level high and low of clean water tank		16.1 liter,Stainless Steel 304 material water level high and low of clean water tank				Prime 30/S6 specification(DHF 3-1-1)
7.05	57	fill water	automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor		automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor		automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor		automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor		automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor		automatic fill water for clean water tank, water pressure 2-7 bar automatic fill water for steam generator by fill water motor				Prime 30/S6 specification(DHF 3-1-1)
7.06	58	Drain	quick drain port		quick drain port		quick drain port		quick drain port		quick drain port		quick drain port				Prime 30/S6 specification(DHF 3-1-1)
8		Air removes / dry / insulation															
8.01	59	gravity (no vacuum pump type)	gravity pulse		gravity pulse		gravity pulse		gravity pulse		gravity pulse		gravity pulse				Prime 30/S6 Operation Manual(DHF 4-2-2.a)
8.02	60	Pre Vacuum	NA		vacuum pump spec. Max Flow XX LPM, MAX Vacuum -0.XX bar Pulse 2 times Max. vacuum value -0.9 bar		vacuum pump spec. Max Flow XX LPM, MAX Vacuum -0.XX bar Pulse 2 times Max. vacuum value -0.9 bar		NA		vacuum pump spec. Max Flow XX LPM, MAX Vacuum -0.XX bar Pulse 2 times Max. vacuum value -0.9 bar		vacuum pump spec. Max Flow XX LPM, MAX Vacuum -0.XX bar Pulse 2 times Max. vacuum value -0.9 bar				Prime 30/S6 specification(DHF 3-1-1)
8.03	61	Post Vacuum	NA		Yes		Yes		NA		Yes		Yes				Prime 30/S6 specification(DHF 3-1-1)
8.04	62	Heater	Bend-Heater 900W		Bend-Heater 900W		Bend-Heater 900W		Bend-Heater 900W		Bend-Heater 900W		Bend-Heater 900W				Prime 30/S6 specification(DHF 4-2-1)
8.05	63	Insulation	keep warm 0h~24h Temperature 45-60°C		keep warm 0h~24h Temperature 45-60°C		keep warm 0h~24h Temperature 45-60°C		keep warm 0h~24h Temperature 45-60°C		keep warm 0h~24h Temperature 45-60°C		keep warm 0h~24h Temperature 45-60°C				Prime 30/S6 Operation Manual(DHF 4-2-2.a)
9		AIR FILTER															
9.01	64	Vacuum remove air filter	HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%				Prime 30/S6 specification(DHF 4-2-1)
9.02	65	Exhaust Filter	NA		NA		filter chamber material SUS304 sterilization 150 cycle times @ 134C temperature sensor K type		NA		NA		filter chamber material SUS304 sterilization 150 cycle times @ 134C temperature sensor K type				Prime 30/S6 specification(DHF 4-2-1)
9.03	66	High pressure cooling air filter	NA		NA		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%		NA		NA		HEPA air filter 0.03um EN 13060, BFE and VFE filtration efficiency as high as 99.999%				Prime 30/S6 specification(DHF 4-2-1)
10		COOLING SYSTEM															
10.01	67	Exhaust	Adjustable Level 0 - 10		Adjustable Level 0 - 10		Adjustable Level 0 - 10		Adjustable Level 0 - 10		Adjustable Level 0 - 10		Adjustable Level 0 - 10				Prime 30/S6 specification(DHF 4-2-1)
10.02	68	Fan cooling	DC 24V x 1 pcs		DC 24V x 3 pcs		DC 24V x 3 pcs		DC 24V x 1 pcs		DC 24V x 1 pcs		DC 24V x 3 pcs				Prime 30/S6 specification(DHF 4-2-1)
10.03	69	Condenser	NA		NA		NA		NA		NA		NA				
10.04	70	High pressure cooling	NA		NA		pressure adjustment valve (with air-water separation), @ 2.0 bar piping cooling rate		NA		NA		pressure adjustment valve (with air-water separation), @ 2.0 bar piping cooling rate				Prime 30/S6 specification(DHF 4-2-1)
11		Data Recorder & Traceable System															

Parameter for Prime 30L, 56L										
No.		Item	Prime 30L Economic	Prime 30L Standard	Prime 30L Advance	Prime 56L Economic	Prime 56L Standard	Prime 56L Advance	Remark / Reference	
11.01	71	Printer	NA	Integrated thermal printer	Integrated thermal printer	NA	Integrated thermal printer	Integrated thermal printer	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
11.02	72	Printer Paper	NA	Thermal paper x 5 pcs	Thermal paper x 5 pcs	NA	Thermal paper x 5 pcs	Thermal paper x 5 pcs	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
11.03	73	Memory	USB memory max. 64GB	USB memory max. 64GB	USB memory max. 64GB	USB memory max. 64GB	USB memory max. 64GB	USB memory max. 64GB	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
11.04	74	Recorder reading software	FileChecker	FileChecker	FileChecker	FileChecker	FileChecker	FileChecker	<a href="#">Prime 30/56 Operation Manual(DHF 4-2-2.a)</a>	
11.05	75	Traceable mode	Pre-scan and post scan scanner bar code printer (@150℃)	Pre-scan and post scan scanner bar code printer (@150℃)	Pre-scan and post scan scanner bar code printer (@150℃)	Pre-scan and post scan scanner bar code printer (@150℃)	Pre-scan and post scan scanner bar code printer (@150℃)	Pre-scan and post scan scanner bar code printer (@150℃)	<a href="#">Prime 30/56 specification(DHF 3-1-1)</a>	
11.06	76	IoT server	IoT connect. by wifi	IoT connect. by wifi or RJ-45	IoT connect. by wifi or RJ-45	IoT connect. by wifi	IoT connect. by wifi or RJ-45	IoT connect. by wifi or RJ-45	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12		Accessories								
12.01	77	basic panel	Standard	Standard	Standard	Standard	Standard	Standard	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.02	78	software	Standard	Standard	Standard	Standard	Standard	Standard	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.03	79	Wire-Mesh Basket	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.04	80	Waste Sterilization Box	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.05	81	Tray Set	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.06	82	Barcode Scanner	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.07	83	Barcode Printer	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.08	84	Water Distiller	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.09	85	RO Water Filter	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.10	86	Wi-Fi Dongle	Optional	Optional	Optional	Optional	Optional	Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
12.11	87	SteriProcess IoT Server	Standard / Optional	Standard / Optional	Standard / Optional	Standard / Optional	Standard / Optional	Standard / Optional	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
13		Documents:								
		<a href="#">User manual</a>								
13.01	88	Catalog with Technical Sheet	Yes	Yes	Yes	Yes	Yes	Yes		
13.02	89	Operation Manual	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 Operation Manual(DHF 4-2-2.a)</a>	
13.03	90	Service Manual	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 Service Manual(DHF 4-2-2.b)</a>	
13.04	91	Warranty of manufacture defect (since	one Year	one Year	one Year	one Year	one Year	one Year		
13.05	92	Spare part (accessories) Manual	Yes	Yes	Yes	Yes	Yes	Yes		
13.06	93	Availability of providing spare parts for	Yes	Yes	Yes	Yes	Yes	Yes		
		Provide product operation training, preventive maintenance, frequently defect solution to the operator	Yes	Yes	Yes	Yes	Yes	Yes		
14		Standards and Directives								
14.01	95	ISO 13485-2016 Quality Management System - Medical Devices	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">ISO 13485 Cer.</a>	
14.02	96	PED 2014/68/EU Pressure Equipment Directive	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">PED Cer.</a>	
14.03	97	ASME BPVC Section VIII-1 Rules for Construction of Pressure Vessels	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 specification(DHF 4-2-1)</a>	
14.04	98	EN/IEC 61010-1 Safety Requirements for measurement control and laboratory use	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 LVD Cer.(DHF 4-8-2)</a>	
14.05	99	EN/IEC 61010-2-040 Safety requirements for sterilizers used to treat medical materials	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 LVD Cer.(DHF 4-8-2)</a>	
14.06	100	EN/IEC 61326-1 EMC	Yes	Yes	Yes	Yes	Yes	Yes	<a href="#">Prime 30/56 EMC Cer.(DHF 4-8-1)</a>	
14.07	101	RoHS Restriction of Hazardous Substances Directive 2011/65/EU	Yes	Yes	Yes	Yes	Yes	Yes	N/A	

Prime 30L Cycle Time				
Program	Loading		Condition	Time
Solid	N/A	121	S : 20 min 、 Dry 20 min	02:27:28
		134	S : 20 min 、 Dry 20 min	01:57:19
Agar	250ml*3	121	S : 30 min 、 EX10 、 Ins.60 、 1H 、 CT80 、 FT95	02:42:30
Dissolution	250ml*3	80	S : 15 min 、 CT80 、 FT95	00:36:30
B&D				02:06:27
Flash	7KG			00:40:41
Latex	7KG			01:30:55
Liquid	250ml*3	121	S : 30min 、 EX10 、 CT80 、 FT95	01:44:29
Liquid	250ml*3	134	S : 15min 、 EX0 、 CT80 、 FT95	01:20:30
Waste Solid	7KG			01:48:31
Waste Mixed	250ml*1+1KG			01:27:05
Waste Liquid				01:35:17
USER	7KG		S : 10min 、 EX0 、 CT80 、 FT95	01:10:54