Autoclawe

SAP Series

SAP are automatic pressure control autoclaves. Its advantage lies in its product line diversity and flexible options to make the right match and perfect fit for your organization.





SAP Series

Features

SAP Series

- n Independent steam generator
- n Cycle progress indicator light
- n Easily accessible emergency stop button
- n Special liquid program
- n Emergency button
- n Cylindrical chamber volume: 100, 143, 196, 350 (L)
- n Square Chamber volume: 110, 260, 454, 680, 848 (L)





SAP Series

SAP Series External View : Loading door

SAP

A Temperature gauge of chamber



B Door handle



C It can change into square door.





SAP Series

SAP Series External View : Unloading door

SAP SERIES WITH DOUBLE DOOR- SAP-400D/450D /500D /600D SAP-S0110D/S0260D/S0454D/S0680D/S0848D

A Door handle



B It can change into square door.








Autoclaw

SAP Series


SAP Series Control Panel : Loading door





SAP Series

SAP Series Control Panel : Unloading door





Certifications

SAP Series

- ISO 13485:2016 Quality Management System Medical Devices
 PED 2014/68/EU Pressure Equipment Directive



SAP Series

SAP Series Specifications

SPEC/Model NO.	SAP-400	SAP-450	SAP-500	SAP-600	SAP-S0110	SAP-S0260	SAP-S0454	SAP-S0680	SAP-S0848	
Chamber Capacity (L)	100	143	196	350	110	260	454	680	848	
Chamber Size	400 × 000	/ ₁ E0 × 000	F00 v 1000	610 v 1200	/.10 v /.10 v cc0	E10 v E10 v 1000	610 v 610 v 1000	610 v 015 v 1220	610 4015 41520	
Diam x Depth (mm)	400 x 800	450 x 900	500 x 1000	610 x 1200	410 x 410 x 660	510 x 510 x 1000	610 x 610 x 1220	610 x 915 x 1220	610 x 915 x 1520	
Overall Dimensions										
WxHxD(mm)	800x1800x1150	800x1800x1250	800x1800x1350	900x1800x1550	950x1800x1350	950x1800x1450	1050x1800x1650	1740x2060x1960	1740x2100x2300	
Voltage (V)				230)/380(Mono/Tri ph	ase)				
Heater (KW)	-	7	9	12	9	12	15 x 2	15 x 2	16 x 2	
Water Fill *		Piped								
Chamber Type		Cylin	drical		Square					
Pressure Option (bar)					0.9-2.1					
Selectable Sterilization					0 60 min					
Timer		0-60 min.								
Selectable Dry Timer					0-60 min.					
Digital Printer					N/A					
6-Point Temp										
Sensor Recorder	0									
1-Point Temp					0					
Sensor Recorder										
Liquid		S								

^{* &}quot;Piped" - uses local water plumbing with auto-fill.



SAP Series

SAP Series (With Double Door) Specifications

SPEC/Model NO.	SAP-400D	SAP-450D	SAP-500D	SAP-600D	SAP-S0110D	SAP-S0260D	SAP-S0454D	SAP-S0680D	SAP-S0848D	
Chamber Capacity (L)	100	143	196	350	110	260	454	680	848	
Chamber Size	/100 v 000	/ ₁ F0 × 000	F00 v 1000	610 v 1200	/.10 v /.10 v 660	F10 v F10 v 1000	610 v 610 v 1000	610 v 015 v 1220	610 v 015 v 1520	
Diam x Depth (mm)	400 x 800	450 x 900	500 x 1000	610 x 1200	410 x 410 x 660	510 x 510 x 1000	610 x 610 x 1220	610 x 915 x 1220	610 x 915 x 1520	
Overall Dimensions										
$W \times H \times D (mm)$	800x1700x1200	800x1700x1300	800x1700x1350	900x1800x1540	850x1710x1110	1025x1795x1480	1170x1850x1800	1740x2073x1810	1740x2073x2130	
Voltage (V)				230	0/380(Mono/Tri ph	ase)				
Heater (KW)		7	9	12	9	12	15 x 2	15 x 2	16 x 2	
Water Fill *	Piped									
Chamber Type		Cylin	drical		Square					
Pressure Option (bar)					0.9-2.1					
Selectable Sterilization					0 60 min					
Timer		0-60 min.								
Selectable Dry Timer					0-60 min.					
Digital Printer					N/A					
6-Point Temp										
Sensor Recorder		0								
1-Point Temp					0					
Sensor Recorder	U									
Liquid		S								

^{* &}quot;Piped" - uses local water plumbing with auto-fill.



SAP Series

SAP Series Specifications

SPEC/Model NO.	SAP-400	SAP-450	SAP-500	SAP-600	SAP-S0110	SAP-S0260	SAP-S0454	SAP-S0680	SAP-S0848
Temperature Indicator					Gauge				
Pressure Indicator *2					Gauge				
Progress Indicator (bar or lights)					S				
Power Indicator					S				
Pressure Auto Lock					0				
Reset					S				
Standby * ³					N/A				
Program Selector					S				
Pre-Vacuum					N/A				
Post Dry Vacuum					0				
Water Pump					0				

Note: "S" - Standard configuration. "0" - Optional configuration.

Other notes:

- Standard Chamber / Jacket / Door / Housing are made of Stainless Steel #304. #316 Stainless steel optionally avaiable for Chamber / Jacket / Door.
- Heater material is Stainless Steel #316.
- Standard Accessories: Tray*1, Glass Viewer Tube * 2
- Optional Accessories: Trolley / Plate Set / Basket
- Sturdy also offers customized sizes up to 2000 L.



^{*&}lt;sup>2</sup> "Gauge" - Includes both Jacket and Chamber pressure gauges.

^{*&}lt;sup>3</sup> "Standby" - Allows the Big Autoclave to maintain jacket temperatures betweeen sterilization sessions.

SAP Series

SAP Series (With Double Door) Specifications

SPEC/Model NO.	SAP-400D	SAP-450D	SAP-500D	SAP-600D	SAP-S0110D	SAP-S0260D	SAP-S0454D	SAP-S0680D	SAP-S0848D
Temperature Indicator					Gauge				
Pressure Indicator *2					Gauge				
Progress Indicator (bar or lights)					S				
Power Indicator					S				
Pressure Auto Lock					0				
Reset					S				
Standby * ³					N/A				
Program Selector					S				
Pre-Vacuum					N/A				
Post Dry Vacuum					0				
Water Pump					0				

Note: "S" - Standard configuration. "0" - Optional configuration.

Other notes:

- Standard Chamber / Jacket / Door / Housing are made of Stainless Steel #304. #316 Stainless steel optionally avaiable for Chamber / Jacket / Door.
- Heater material is Stainless Steel #316.
- Standard Accessories: Tray*1, Glass Viewer Tube * 2
- Optional Accessories: Trolley / Plate Set / Basket
- Sturdy also offers customized sizes up to 2000 L.



^{*2 &}quot;Gauge" - Includes both Jacket and Chamber pressure gauges.

^{*&}lt;sup>3</sup> "Standby" - Allows the Big Autoclave to maintain jacket temperatures betweeen sterilization sessions.

SAP Series

Comparison Table - Floor Stand Series Feature Comparison

Model NO. Specification	HP SAT		SAP	Pointer-Auto	Pointer-Manual				
Control System	Digital Temp. Pl	.C Control	Mechanical Pressure Control Optional: (T) for Digital Temp. Control	Mechanical P	Pressure Control				
Program for selection Wrapped/Unwrapped/Liquid		Yes		A					
Air Removal	Vacuum P	ump		Gravity					
Pre-Vacuum Function	Yes			N/A					
Post(Dry)-Vacuum Function	Yes		Optional (V)		N/A				
Dry function			Yes (Auto)	es (Auto) Yes (Manual)					
Printer & USB storage	Yes								
1-point temp. recorder	N/A	Yes	Optional (X)		/A				
6-point temp. recorder	N/A		Optional (Y)	N/A					
Double Door		Optional ([0)	N/A					
Fill water system*		Auto	٦	Manual Optional (B) for auto	Manual				
(Water pressure 2 bar required)		, , , , , , , , , , , , , , , , , , ,	S	riariaat					
Chamber/Jacket	Gauge / LCD Display		Gauge						
Pressure display	oddyc / Lob bispidy		dayo						
Stand-by Switch	Yes								
Emergency/Reset Switch		Ye	S	N/A					
Add-Water Pump		Option	al (W)	Optional (BW)	N/A				
Water level Protection			Yes Optional (L)						
Pressure door auto-lock	Yes		Optional (K)						

Option Codes (e.g. Optional (V)) listed in the table above are only reference codes used by other documentation.



Optional Codes

SAP Series

Option Codes (e.g. Optional (V)) listed in the table above are only reference codes used by other documentation.

D: DOUBLE DOOR DEVICE

K: PRESSURE DOOR AUTO-LOCK DEVICE

T: DIGITAL DISPLAYED TEMP. CONTROLLER

V: DRY VACUUM PUMP

W: ADD WATER PUMP

X: 1-POINT TEMP. RECORDER

Y: 6-POINT TEMP. RECORDER



SAP Series

Accessories





B Plate Set



C.1 Optional Trolley



D Filter



C.2 Optional Trolley



