SAT Series

SAT are digital temperature control autoclaves with pre & dry vacuum. Designed for Hospital medical, laboratory, biological science, pharmacy industry, reliability testing equipment.





SAT Series

Features

SAT Series

- Fractioned pre vacuum air removal
- Independent steam generator
- Built-in temperature recorder
- Cycle progress indicator light
- Double locking safety system
- Special liquid program
- Stand by function
- Cylindrical chamber volume: 100, 143, 196, 350 (L)
- Square Chamber volume: 110, 260, 454, 680, 848 (L)






SAT Series

SAT Series External View : Loading door

SAT

A Recorder







C It can change into square door.





SAT Series

SAT Series External View : Unloading door

SAT SERIES WITH DOUBLE DOOR- SAT-400D/450D /500D /600D

SAT-S0110D/S0260D/S0454D/S0680D/S0848D












SAT Series


SAT Series Control Panel: Loading door



Autoclawe

SAT Series

SAT Series Control Panel : Unloading door





Certifications

SAT Series

- ISO 13485:2016 Quality Management System Medical Devices
 PED 2014/68/EU Pressure Equipment Directive



SAT Series

SAT Series Specifications

SPEC/Model NO.	SAT-400	SAT-450	SAT-500	SAT-600	SAT-S0110	SAT-S0260	SAT-S0454	SAT-S0680	SAT-S0848		
Chamber Capacity (L)	100	143	196	350	110	260	454	680	848		
Chamber Size											
Diam x Depth (mm)	400x800	450x900	500x1000	610x1200	410x410x660	510x510x1000	610x610x1220	610x915x1220	610x915x1520		
Overall Dimensions											
WxHxD(mm)	800x1800x1150	800x1800x1250	800x1760x1440	900x1800x1550	950x1800x1350	950x1800x1550	1450x1925x1800	1740x2060x1960	1740x2060x2260		
Voltage (V)					230/380(Tri phase)					
Heater (KW)	7 9			12	9	12	15 x	2	16 x 2		
Water Fill *		Piped									
Chamber Type		Cylind	rical			Square					
Temperature		105 170									
Options (°C)		105-136									
Selectable											
Sterilization Timer		0-60 min.									
Selectable Dry Timer					0-60 min.						
Digital Printer					N/A						
6-Point Temp		0									
Sensor Recorder											
1-Point Temp					C						
Sensor Recorder		S									
Liquid		S									

^{* &}quot;Piped" - uses local water plumbing with auto-fill.



SAT Series

SAT Series (With Double Door) Specifications

SPEC/Model NO.	SAT-400D	SAT-450D	SAT-500D	SAT-600D	SAT-S0110D	SAT-S0260D	SAT-S0454D	SAT-S0680D	SAT-S0848D		
Chamber Capacity (L)	100	143	196	350	110	260	454	680	848		
Chamber Size											
Diam x Depth (mm)	400x800	450x900	500x1000	610x1200	410x410x660	510x510x1000	610x610x1220	610x915x1220	610x915x1520		
Overall Dimensions											
WxHxD(mm)	800x1700x1300	800x1700x1400	800x1700x1450	900x1800x1640	850x1710x1110	1025x1795x1480	1170x1850x1800	1740x2073x1810	1740x2073x2130		
Voltage (V)					230/380(Tri phase	J					
Heater (KW)	7		9	12	9	12	15 x	2	16 x 2		
Water Fill *					Piped						
Chamber Type		Cylind	rical		Square						
Temperature		105 170									
Options (°C)	105-136										
Selectable											
Sterilization Timer		0-60 min.									
Selectable Dry Timer					0-60 min.						
Digital Printer					N/A						
6-Point Temp		0									
Sensor Recorder											
1-Point Temp					C						
Sensor Recorder		S									
Liquid		S									

^{* &}quot;Piped" - uses local water plumbing with auto-fill.



SAT Series

SAT Series Specifications

SPEC/Model NO.	SAT-400	SAT-450	SAT-500	SAT-600	SAT-S0110	SAT-S0260	SAT-S0454	SAT-S0680	SAT-S0848			
Temperature Indicator					LED display							
Pressure Indicator *2		Gauge										
Progress Indicator					C							
(bar or lights)					S							
Power Indicator		S										
Pressure Auto Lock					S							
Reset					S							
Standby * ³					S							
Program Selector					S							
Pre-Vacuum					S							
Post Dry Vacuum					S							
Water Pump					0							

Note: "S" - Standard configuration. "0" - Optional configuration.

Other notes:

- Standard Chamber / Jacket / Door / Housing are made of Stainless Steel #304. #316 Stainless steel optionally avaiable for Chamber / Jacket / Door.
- Heater material is Stainless Steel #316.
- Standard Accessories: Tray*1, Glass Viewer Tube * 2
- Optional Accessories: Trolley / Plate Set / Basket
- Sturdy also offers customized sizes up to 2000 L.



^{*&}lt;sup>2</sup> "Gauge" - Includes both Jacket and Chamber pressure gauges.

^{*&}lt;sup>3</sup> "Standby" - Allows the Big Autoclave to maintain jacket temperatures betweeen sterilization sessions.

SAT Series

SAT Series (With Double Door) Specifications

SPEC/Model NO.	SAT-400D	SAT-450D	SAT-500D	SAT-600D	SAT-S0110D	SAT-S0260D	SAT-S0454D	SAT-S0680D	SAT-S0848D		
Temperature Indicator					LED display						
Pressure Indicator *2		Gauge									
Progress Indicator (bar or lights)					S						
Power Indicator					S						
Pressure Auto Lock					S						
Reset					S						
Standby *3					S						
Program Selector					S						
Pre-Vacuum					S						
Post Dry Vacuum					S						
Water Pump					0						

Note: "S" - Standard configuration. "0" - Optional configuration.

Other notes:

- Standard Chamber / Jacket / Door / Housing are made of Stainless Steel #304. #316 Stainless steel optionally avaiable for Chamber / Jacket / Door.
- Heater material is Stainless Steel #316.
- Standard Accessories: Tray*1, Glass Viewer Tube * 2
- Optional Accessories: Trolley / Plate Set / Basket
- Sturdy also offers customized sizes up to 2000 L.



^{*2 &}quot;Gauge" - Includes both Jacket and Chamber pressure gauges.

^{*&}lt;sup>3</sup> "Standby" - Allows the Big Autoclave to maintain jacket temperatures betweeen sterilization sessions.

SAT Series

Comparison Table - Floor Stand Series Feature Comparison

Model NO. Specification	HP SAT		SAP	Pointer-Auto	Pointer-Manual		
Control System	Digital Temp. PL	C Control	Mechanical Pressure Control Optional: (T) for Digital Temp. Control	ressure Control			
Program for selection Wrapped/Unwrapped/Liquid		Yes		А			
Air Removal	Vacuum P	ump		Gravity			
Pre-Vacuum Function	Yes			N/A			
Post(Dry)-Vacuum Function	Yes		Optional (V)	Optional (V)			
Dry function			Yes (Auto)				
Printer & USB storage	Yes N/A						
1-point temp. recorder	N/A	Yes	Optional (X)	A			
6-point temp. recorder	N/A		Optional (Y)	N/A			
Double Door		Optional ([
Fill water system* (Water pressure 2 bar required)		Auto		Manual			
Chamber/Jacket Pressure display	Gauge / LCD Display		Gauge	Gauge			
Stand-by Switch	Yes						
Emergency/Reset Switch		Ye	S	N/A			
Add-Water Pump		Optiona	al (W)	Optional (BW)	N/A		
Water level Protection			Yes Optional (L)				
Pressure door auto-lock	Yes		Optional (K)				

Option Codes (e.g. Optional (V)) listed in the table above are only reference codes used by other documentation.



Optional Codes

Option Codes (e.g. Optional (V)) listed in the table above are only reference codes used by other documentation.

D: DOUBLE DOOR DEVICE W: ADD WATER PUMP

Y: 6-POINT TEMP. RECORDER

SAT Series



SAT Series

Accessories





B Plate Set



C.1 Optional Trolley



D Filter



C.2 Optional Trolley

